

THERMCELL MEDIKA

This device combines three methods: Thermolysis, Electrolysis and Blend. ThermCell Medika is a device that perfectly cope with cracked capillaries on the face, remove the single hairs and small skin lesions. If you want to learn more about the innovative method Blend, click below to see the details.

Application

Closing dilated capillares
Deleting single hairs
Removal of minor of skin lesions

TREATMENT PROGRAMS

THERMOLYSIS

Treatment using alternating current for high-frequency and low intensity

ELECTROLYSIS

Treatment using direct current (galvanic)

BLEND

Method using both activity galvanic current and alternating current

BENEFITS

High efficiency
Immediate effects after one treatment
The precision of the device and safety
Full control of the pulses intensity
Disposable needles Balette type

Minimal discomfort - anesthesia is not required
Readable touch screen for easy operation with thermolysis
Intuitive device menu in your native language
The device works only on the spot treatment without damaging surrounding tissue
Elegant black and white design made from Corian®
Modern design
Operating mode continuous or pulsed



