

## **THERMCELL MEDIKA**

This device combines three methods: Thermolysis, Electrolysis and Blend. ThermCell Medika is a device that perfectly cope with cracked capillaries on the face, remove the single hairs and small skin lesions. If you want to learn more about the innovative method Blend, click below to see the details.

### Application

Closing dilated capillares  
Deleting single hairs  
Removal of minor of skin lesions

### TREATMENT PROGRAMS

#### **THERMOLYSIS**

Treatment using alternating current for high-frequency and low intensity

#### **ELECTROLYSIS**

Treatment using direct current (galvanic)

#### **BLEND**

Method using both activity galvanic current and alternating current

### BENEFITS

High efficiency  
Immediate effects after one treatment  
The precision of the device and safety  
Full control of the pulses intensity  
Disposable needles Balette type

Minimal discomfort - anesthesia is not required  
Readable touch screen for easy operation with thermolysis  
Intuitive device menu in your native language  
The device works only on the spot treatment without damaging surrounding tissue  
Elegant black and white design made from Corian®  
Modern design  
Operating mode continuous or pulsed



